



**MEDIHONEY®**  
Wound and Burn Dressing

## Advanced Wound Care



**Day 0**  
DFU measuring 3.0 x 3.0 cm with 100% necrotic tissue.



**Day 2**  
2 days after MediHoney was applied, wound is still 3.0 x 3.0 cm, but now has yellow stringy slough and granulation tissue.



**Week 3**  
After 2 more applications of MediHoney, the wound measured 2.5 x 2.5 cm and had 20% yellow stringy slough and 80% granulation tissue.

## MediHoney® Dressings

Aids and supports autolytic debridement and provides a moist wound healing environment in acute and chronic wounds and burns<sup>1-4</sup>

MediHoney® Dressings, with Active *Leptospermum* Honey, aids and supports autolytic debridement and a moist wound healing environment in acute and chronic wounds and burns.<sup>1-4</sup> This is done through two properties: High Osmolarity (Fig. 1) and Low pH (Fig. 2).

### Features & Benefits:

- MediHoney dressings address many factors that delay healing, helps to promote a moist wound environment that aids and supports autolytic debridement
- Promotes a moisture-balanced wound environment conducive to wound healing in multiple etiologies<sup>5</sup>
- Osmotic potential draws fluid through the wound, to the surface, helping to liquefy non-viable tissue
- MediHoney has a low pH of 3.5-4.5, lowering pH has been associated with wound healing benefits

MediHoney – Committed to bringing you quality, versatility and innovative solutions for your patients' care throughout the continuum.




Fig. 1




Fig. 2

	Reference	Description	Packaging Unit/Case	HCPCS
<b>Gel</b>	31805	0.5 oz tube	10/box, 4 boxes/case	A4649
	31815	1.5 oz tube	1/box, 12 boxes/case	A4649
	31840	14 oz tub	1/jar, 6 tubs/case	-----
<b>Paste</b>	31505	0.5 oz tube	10/box, 4 boxes/case	A4649
	31515	1.5 oz tube	1/box, 12 boxes/case	A4649
	31535	3.5 oz tube	1/box, 12 boxes/case	A4649
<b>Hydrogel Sheet Non-Adhesive</b>	31620	2.4 in x 2.4 in	10/box, 5 boxes/case	<b>A6242</b>
	31640	4.3 in x 4.3 in	10/box, 5 boxes/case	<b>A6243</b>
<b>Adhesive</b>	31720	2.8 in x 2.8 in (4.3 in x 4.3 in with adhesive border)	10/box, 5 boxes/case	<b>A6245</b>
	31740	4.5 in x 4.5 in (6 in x 6 in with adhesive border)	10/box, 5 boxes/case	<b>A6246</b>
<b>HCS Surgical</b>	31738	1.75 in x 6.5 in (3 in x 8 in with adhesive border)	10/box, 5 boxes/case	A4649
<b>Fenestrated (Non-Adhesive)</b>	31618	1.8 in x 1.8 in	10/box, 5 boxes/case	A4649
	31622	2.4 in x 2.4 in	10/box, 5 boxes/case	A4649
	31644	4.33 in x 4.33 in	10/box, 5 boxes/case	A4649
<b>Non-Adhesive</b>	31612	8 in x 12 in	2/box, 5 boxes/case	A4649
	31722	2.8 in x 2.8 in (4.3 in x 4.3 in with adhesive border)	10/box, 5 boxes/case	A4649
	31744	4.5 in x 4.5 in (6 in x 6 in with adhesive border)	10/box, 5 boxes/case	A4649
<b>Adhesive</b>	31722	2.8 in x 2.8 in (4.3 in x 4.3 in with adhesive border)	10/box, 5 boxes/case	A4649
	31744	4.5 in x 4.5 in (6 in x 6 in with adhesive border)	10/box, 5 boxes/case	A4649
	31012	0.75 in x 12 in	5/box, 4 boxes/case	A4649
<b>Calcium Alginate</b>	31012	0.75 in x 12 in	5/box, 4 boxes/case	A4649
	31022	2 in x 2 in	10/box, 10 boxes/case	A4649
	31045	4 in x 5 in	10/box, 5 boxes/case	A4649



Gel



Paste



Hydrogel Sheet  
(Non-Adhesive)



Hydrogel Sheet  
(Adhesive)



HCS (Surgical)



HCS - Fenestrated  
(Non-Adhesive)



HCS (Non-Adhesive)



HCS (Adhesive)



Calcium Alginate

References: 1. Regulski, M. A novel wound care dressing for chronic leg ulcerations. Podiatry Management, 2008. November/December: p. 235-246 2. Robson, V., Dodd, S and Thomas, S. Standardized antibacterial honey (MediHoney®) with standard therapy in wound care: randomized clinical trial. Journal of Advanced Nursing, 2009; p. 565-575. 3. Bateman S, Graham T (2007) The Use of MediHoney® Wound Gel on surgical wounds post-CABG. WOUNDS UK; Vol 3(3). 76 – 83. 4. Cadogan, J. (2008) The use of honey to treat an ulcer on the heel of a person with diabetes. The Diabetic Foot Journal;11. (1): 43-45. 5. in house data.

Integra LifeSciences Corporation intends to use reasonable efforts to provide accurate coding information, but this information should not be construed as providing clinical advice, dictating reimbursement policy or substituting for the judgment of a practitioner. It is always the Provider's responsibility to determine and submit appropriate codes, charges and modifiers for services that are rendered. Integra LifeSciences Corporation assumes no responsibilities or liabilities for the timeliness, accuracy and completeness of the information contained herein. Since reimbursement laws, regulations and payor policies change frequently, it is recommended that providers consult with their payors, coding specialists and/or legal counsel regarding coverage, coding and payment issues.

**Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.**

- Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- Warning: Applicable laws restrict these products to sale by or on the order of a physician.
- Consult product labels and inserts for any indications, contraindications, hazards, warnings, precautions, and instructions for use.

**For more information or to place an order, please contact:**

**United States**

USA 800-654-2873 ▪ 888-980-7742 fax

International +1 609-936-5400 ▪ +1 609-750-4259 fax

[integralife.com](http://integralife.com)

