



MEDIHONEY®
Wound and Burn Dressing

PREPARE


Advanced Wound Care



Gel



100% Paste



HCS
(Hydrogel Colloidal Sheet)



Calcium Alginate

MediHoney® Dressings


Aids and supports
autolytic debridement
and a moist wound
healing environment
in acute and chronic
wounds and burns¹⁻⁴

MediHoney® Dressings, made with Active *Leptospermum* Honey, offer versatility and demonstrated effectiveness when managing your patients' superficial to deep partial thickness burns⁵⁻⁹

Features & Benefits:

- Promotes a moisture-balanced wound environment conducive to wound healing in 1st and 2nd degree burns¹⁰
- Osmotic potential draws fluid through the wound, to the surface, helping to liquefy non-viable tissue
- MediHoney has a low pH of 3.5-4.5, lowering pH has been associated with wound healing benefits
- Supports the moist wound healing of deep partial thickness burns when you prefer not to graft*

MediHoney – Committed to bringing you quality, versatility and innovative solutions for your patients' care throughout the continuum.



*Cleared indications for use: 1st and 2nd degree burns.

INTEGRA®
LIMIT UNCERTAINTY

	Reference	Description	Packaging Unit/Case	HCPCS
Gel	31805	0.5 oz tube	10/box, 4 boxes/case	A4649
	31815	1.5 oz tube	1/box, 12 boxes/case	A4649
	31840	14 oz tub	1/jar, 6 tubs/case	-----
Paste	31505	0.5 oz tube	10/box, 4 boxes/case	A4649
	31515	1.5 oz tube	1/box, 12 boxes/case	A4649
	31535	3.5 oz tube	1/box, 12 boxes/case	A4649
Hydrogel Sheet Non-Adhesive	31620	2.4 in x 2.4 in	10/box, 5 boxes/case	A6242
	31640	4.3 in x 4.3 in	10/box, 5 boxes/case	A6243
Adhesive	31720	2.8 in x 2.8 in (4.3 in x 4.3 in with adhesive border)	10/box, 5 boxes/case	A6245
	31740	4.5 in x 4.5 in (6 in x 6 in with adhesive border)	10/box, 5 boxes/case	A6246
HCS Surgical	31738	1.75 in x 6.5 in (3 in x 8 in with adhesive border)	10/box, 5 boxes/case	A4649
Fenestrated (Non-Adhesive)	31618	1.8 in x 1.8 in	10/box, 5 boxes/case	A4649
	31622	2.4 in x 2.4 in	10/box, 5 boxes/case	A4649
	31644	4.33 in x 4.33 in	10/box, 5 boxes/case	A4649
Non-Adhesive	31612	8 in x 12 in	2/box, 5 boxes/case	A4649
	31722	2.8 in x 2.8 in (4.3 in x 4.3 in with adhesive border)	10/box, 5 boxes/case	A4649
	31744	4.5 in x 4.5 in (6 in x 6 in with adhesive border)	10/box, 5 boxes/case	A4649
Adhesive	31722	2.8 in x 2.8 in (4.3 in x 4.3 in with adhesive border)	10/box, 5 boxes/case	A4649
	31744	4.5 in x 4.5 in (6 in x 6 in with adhesive border)	10/box, 5 boxes/case	A4649
	31012	0.75 in x 12 in	5/box, 4 boxes/case	A4649
Calcium Alginate	31012	0.75 in x 12 in	5/box, 4 boxes/case	A4649
	31022	2 in x 2 in	10/box, 10 boxes/case	A4649
	31045	4 in x 5 in	10/box, 5 boxes/case	A4649



Gel



Paste



Hydrogel Sheet
(Non-Adhesive)



Hydrogel Sheet
(Adhesive)



HCS (Surgical)



HCS - Fenestrated
(Non-Adhesive)



HCS (Non-Adhesive)



HCS (Adhesive)



Calcium Alginate

References: 1. Regulski, M. A novel wound care dressing for chronic leg ulcerations. *Podiatry Management*, 2008. November/December: p. 235-246. 2. Robson, V., Dodd, S and Thomas, S. Standardized antibacterial honey (MediHoney®) with standard therapy in wound care: randomized clinical trial. *Journal of Advanced Nursing*, 2009; p. 565-575. 3. Bateman S, Graham T (2007) The Use of MediHoney® Wound Gel on surgical wounds post-CABG. *WOUNDS UK*; Vol 3(3). 76 – 83. 4. Cadogan, J. (2008) The use of honey to treat an ulcer on the heel of a person with diabetes. *The Diabetic Foot Journal*;11. (1): 43-45. 5. Wijesinghe M, Weatherall M, Perrin K, Beasley R. Honey in the treatment of burns: a systematic review and meta-analysis of its efficacy. *The New Zealand Medical Journal*. 2009; 122(5): 47-60. 6. Becker B. If at first you don't succeed...try HONEY for wound burn resolution. *Approaches in Critical Care*. 2012; 7 (2): 2-4. 7. Subrahmanyam M. Honey Dressing for Burns-An Appraisal. *Annals of Burns and Fire Disaster*. 1996; 9 (1). 8. Puljak L, et al. Commercial dressings for burns versus sweet ancient remedy. *Intern Emerg Med*. 2009; (4):53-56. 9. Subrahmanyam M. A prospective randomised clinical and histological study of superficial burn wound healing with honey and silver sulfadiazine. *Burns*. 1998; 24: 157-161. 10. In house data.

Integra LifeSciences Corporation intends to use reasonable efforts to provide accurate coding information, but this information should not be construed as providing clinical advice, dictating reimbursement policy or substituting for the judgment of a practitioner. It is always the Provider's responsibility to determine and submit appropriate codes, charges and modifiers for services that are rendered. Integra LifeSciences Corporation assumes no responsibilities or liabilities for the timeliness, accuracy and completeness of the information contained herein. Since reimbursement laws, regulations and payor policies change frequently, it is recommended that providers consult with their payors, coding specialists and/or legal counsel regarding coverage, coding and payment issues.

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

- Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- Warning: Applicable laws restrict these products to sale by or on the order of a physician.
- Consult product labels and inserts for any indication, contraindications, hazards, warnings, precautions, and instructions for use.

For more information or to place an order, please contact:

United States, Canada, Asia, Pacific, Latin America

USA 800-654-2873 ■ 888-980-7742 fax

International +1 609-936-5400 ■ +1 609-750-4259 fax

integralife.com

